

# Record an ECG

Version: 1 | Last Modified on 06/05/2026 9:50 am AEST

Zedmed on-prem supports WelchAllyn and OfficeMedic ECGs machines and connects to them using Cardioview from a patient's clinical record. Zedmed Cloud does not support ECG integrations.

To record an ECG

1. From Zedmed Clinica, select on the **Open Patient** Icon.
2. Type in the first few letters of the patient surname, comma, the patient initial.
3. Select **Find**.
4. Select the correct patient, click on **Select**.
5. Select **Perform ECG** icon.



6. This launches the third-party ECG software - **Cardioview**.
  7. The patient's name should be populated.
  8. Perform the ECG as per usual.
  9. Select **Save**.
- A window asking for a description will appear.
10. Enter the title you wish the image to be saved as in the patient record.
  11. Select **OK**.
  12. The saved ECG will be listed on the **Images/ECG** tab in the Summary Views area.
  13. To view the ECG, right-click on it and select **View Image**.
-